

Catheters

Word Bartholin Catheter Set

This latex-free Bartholin Catheter Set contains a silicone inflatable catheter, scalpel, syringe, and needle, for the treatment of Bartholin gland abscesses.

SKU	CPT Code	HCPT Code	
WC-400MM	56420	A4550	6units/pack

Uterine Manipulator Injector

Ideally designed for Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Minilaparotomy, Salpingoplasty, and Fertility Exams.

1. Pear shaped balloon for optimum hold and better seal at the internal OS
2. Easy insertion by reducing the size of cannula
3. Flexible sleeve keeps luer from cracking
4. No need for another instrument during the procedure

SKU	Description	O.D.	
UMI 4.5	Uterine Manipulator/ Injector	4.5mm	12units/box
UMI 4.5K	Uterine Manipulator/ Injector with clamp	4.5mm	12units/box

Uterine Injector

1. Sterile, disposable, ready to use
2. Curved Shaft (uterine sound shaped) for use in anteverted or retroverted uteri
3. Snap-on cervical stop custom set for every patient
4. Snap-on cervical stop V-shaped design for easy view of external cervical OS
5. Patient ambulatory once instrument properly inserted
6. Durable no-crack luer fitting on flexible collar for convenient injection
7. No spring or weights added

SKU	Description	Tube O.D.	
UI 4.0	Uterine Injector 4.0mm	4.0 mm	12units/box

HSG Stylet

SKU	Description	Length	
BA001	H/S Catheter Stylet	420mm	10units/box

HSG Catheter Sets Latex Free / with balloon

Bioteque's latex-free HSG Catheter Sets, available in 5F & 7F with a silicone balloon, are useful in evaluating the causes of abnormal bleeding, menstrual disorders, recurrent miscarriage, and unexplained infertility.

SKU	Description	
BA05F	5 French Catheters	10 units/ box
BA07F	7 French Catheters	10 units/ box
BA051F	5 French SureFlex Catheters	10 units/ box
BA071F	7 French SureFlex Catheters	10 units/ box

Bellota™ Sonohysterography Catheter

The balloonless Sono Catheter is a sterile disposable 5F catheter that consists of an insertion sheath and cervical cone. The adjustable cervical cone positioner helps with catheter placement while reducing leakage of saline at the external cervical OS.

SKU	Description	Length	
UIC5F	SonoCatheter w/o cervical cone 5Fr	327.5mm	25units/box
UIC5FA	SonoCatheter with cervical cone 5Fr	276.5mm	25units/box
UIC5FAX	SonoCatheter with cervical cone 5Fr	259.5mm	25units/box
UIC5FS	SonoCatheter/Uterine Introducing 5Fr	350.0mm	25units/box

CPT Code : 76831

