


Dilators



Silicone Vaginal Dilator Sets



Each set includes 4 graduated sizes of dilators used to enlarge the vaginal for examination, surgical procedures, or relief of vaginismus. Single patient use only.

Large Set

Family	Length	OD Small	OD Large
DT-A/XL	150.0	30.0	38.0
DT-A/L	138.0	25.5	32.0
DT-A/M	128.0	22.5	28.0
DT-A/S	120.0	19.5	24.0

Medium Set


Family	Length	OD Small	OD Large
DT-B/XL	95.0	26.0	33.0
DT-B/L	86.0	22.5	28.0
DT-B/M	78.0	19.5	23.0
DT-B/S	75.0	12.5	17.0

Small Set

Family	Length	OD Small	OD Large
DT-C/XL	75.0	13.5	18.0
DT-C/L	65.0	12.0	16.0
DT-C/M	60.0	10.5	14.0
DT-C/S	50.0	9.5	13.0

Unit : mm

Bioteque Vaginal Dilator: CPT Code: 99070 Procedure Code: 57400



OS Finder Sets 25/Box

Bioteque offers disposable (blue) OS Finder set. The OS Finder set include an OS Finder, Canal Locator, and Fundus Dilator. The disposable OS Finders can be ordered separately. These sets are useful in locating and dilating the external and internal OS as well as the cervical canal.

SKU	Description	Length	Taper	O.D.
OS265BA-1	Disposable OSFinder	208mm	23mm	Ø 1.75 - 3.8mm
OS265BA-2	Disposable Canal Locator	208mm	45mm	Ø 1.75 - 3.8mm
OS265BA-3	Disposable Fundus Dilator	208mm	85mm	Ø 1.75 - 3.8mm



Uterine Sound 25/Box

Bioteque's latex-free Uterine Sounds are malleable, and disposable. Used for dilating and measuring the cervix. Graduated in centimeters with markings at 3cm, 5cm, 7cm, 9cm, and 11cm.

SKU	Description	Length	O.D.
BA100	Disposable Uterine Sound	248.5mm	3mm